

Blood Agar for the Isolation of *M. tuberculosis* in Comparison with LJ Medium and BACTEC MGIT 960

SWETHA KAMBLE¹, P SHASHIKALA REDDY², CH NAVANEETHA³

ABSTRACT

Introduction: The evolution of laboratory diagnostic methods of Tuberculosis (TB) over the past few years are designed to achieve more rapid, less expensive, and accurate results. However, acid-fast staining and culture for Mycobacteria remains the core of any diagnostic algorithm for TB.

Aim: To assess 7% Sheep Blood Agar (BA) for the primary isolation of *Mycobacterium tuberculosis* in comparison with the Lowenstein-Jensen (LJ) medium and BACTEC Mycobacterial Growth Inhibitor Tube (MGIT) 960.

Materials and Methods: This is a prospective study and was carried over a period of seven months in State TB Demonstration and Training Centre Intermediate Reference Laboratory (STDC-IRL), Erragadda Hyderabad, India from March 2014 to September 2014. About 100 clinically suspected TB cases (95 pulmonary and 5 extra pulmonary) were selected and respective samples were subjected to Ziehl-Neelsen (ZN) staining. Twenty one samples were positive and 79 were negative for Acid-Fast Bacilli (AFB). The samples were processed and inoculated on BA slants, LJ medium and BACTEC MGIT 960 (Becton,

Dickinson and Company) bottles. The cultures were confirmed by ZN staining and speciated by biochemical reactions. Further, validation of growth on BA was identified by Genotype MTBDR plus (Hains Lifescience, Nehren, Germany) to confirm the growth of *M. tuberculosis*.

Results: Approximately 100 samples inoculated, 40 {38 *M. tuberculosis* and 2 Non *tuberculous mycobacteria* (NTM)} were isolated on MGIT 960, 37 (34 *M. tuberculosis* and 3 NTM) on BA and 32 (30 *M. tuberculosis* and 2 NTM) on LJ medium. The average isolation time of BA LJ and BACTEC MGIT 960 TB was 13.2 days, 26 days and 11.8 days. Nearly 37 blood culture positive samples tested by Genotype MTBDR plus, 34 (91.89%) were confirmed as *Mycobacterium Tuberculosis* Complex (MTBC). The contamination rate was 2% on BA, 6% on LJ and 12% on BACTEC MGIT 960.

Conclusion: Primary isolation of *M. tuberculosis* was achieved 8 to 18 days after inoculation of clinical samples on BA. It requires no specific equipment and therefore, BA can be prepared, utilised and also have the advantage of shorter turnaround time.

Keywords: Extrapulmonary, Lowenstein jensen medium, Mycobacterial growth inhibitor tube

INTRODUCTION

According to Revised National Tuberculosis Control Programme (RNTCP) annual report, one-fourth of the global TB burden is seen in India. In 2015, 28 lacs cases were estimated and 4.8 lacs people died from TB in India [1].

Various new modalities like Cartridge Based-Nucleic Acid Amplification Techniques (CB-NAAT), Microscopically Observed Drug Susceptibility (MODS), Loop Mediated Isothermal Amplification (LAMP), Line Probe Assay (LPA), Xpert MTB/RIF are employed in the diagnosis of TB [2]. But they are costly; consumables are not easily available and require technical expertise which can limit the utilisation of automated methods. Also the introduction of automated methods does not mean to replace the need for culture as mycobacterial culture is still needed for [1]. Surveillance of drug resistance as part of control programme performance [2] follow up of patients who fail on treatment regimens to check for drug resistance for drug susceptibility and identification tests which include both genotypic and phenotypic tests [3].

In developing countries like India, LJ medium is the most commonly used medium for the culture of *M. tuberculosis*. Growth of *M. tuberculosis* is slow, appear in two to six weeks and negative culture report cannot be given before eight weeks [4]. However few research papers published such as Drancourt M and Raoult D, Satti L et al., Coban AY et al., suggest that BA can also be used for the isolation of *M. tuberculosis* [5-7].

In the present study, an attempt has been made to assess the utility of BA for isolation of *M. tuberculosis* with LJ medium and BACTEC MGIT 960. To make it a selective medium, antibiotics and antifungal drugs were added. Further identification of the isolates recovered from BA were confirmed by Genotype MTBDR plus assay (Hains Lifescience, Nehren, Germany) targeting 23rRNA, which was available at the STDC-IRL.

If this study succeeds in proving the *M. tuberculosis* culture on BA as a possible modality in TB diagnostics, then it will be a path breaking finding in this field.

MATERIALS AND METHODS

In this prospective study, samples were collected from the initial screened 100 clinically suspected TB cases attending Designated Microscopy Centre (DMC) of Government. General and Chest Hospital based on certain inclusion and exclusion criteria as per RNTCP guidelines [3] and were further processed in the STDC-IRL.

Inclusion Criteria

The study group includes patients with more than any two symptoms out of the following: 1. Fever and cough with expectoration >2 weeks; 2. Abnormal findings in chest radiograph; 3. Gradual weight loss; 4. HIV status (positive/negative); 5. Contacts of smear positive TB patients having cough of any duration. A patient with extra

pulmonary TB may have general symptoms like weight loss, fever with evening rise and night sweats. Other symptoms depend on the organ affected [3].

Exclusion Criteria

1. Cases already on Anti Tuberculous Treatment (ATT);
2. Cases already confirmed as TB by other microbiological methods [3].

It included 95 sputum, 3 pleural fluid, 1 pus and 1 lymph node sample. Two sputum samples were collected from each patient, on the spot and early morning sputum sample. The extra pulmonary samples were collected by the treating physician under strict aseptic conditions. The informed consent of the patients was taken before the collection of the specimens.

The samples were then decontaminated and concentrated by using N-acetylcysteine Sodium Hydroxide method [3] for sputum samples and sulfuric acid [8] for extra pulmonary samples. If any delay in the processing, the samples were stored in the refrigerator at 4°C. A microscopic examination (ZN staining) followed by the parallel inoculation on 7% Sheep BA, LJ and BACTEC MGIT 960 (Becton, Dickinson and Company) was done. LJ medium were prepared as per RNTCP guidelines [3]. BA slants were prepared according to Mathur ML et al., instead of using petri dishes; McCartney bottles were used for BA to avoid desiccation [4]. About 50 µL of concentrated sample was inoculated on: (i) two plain LJ slants and one LJ slant with Para Nitro Benzoic (PNB) acid (500 µg/mL); (ii) two BA slants and one BA slant with PNB and 0.5 mL of the inoculum was added to MGIT tubes in duplicate. Growth obtained on solid medium and MGIT tubes were further tested for the presence of AFB by ZN staining and four biochemical tests, i.e., Niacin, Nitrate, Heat resistant catalase test and PNB test to differentiate the growth into *M. tuberculosis* and NTM. Standard controls *Mycobacterium smegmatis* (ATCC700084) and *M. tuberculosis* H37Rv (ATCC19977) were used in the study.

An isolate was designated to be that of *M. tuberculosis*, if it grew slowly (taking more than 7 days), showing no growth on medium containing PNB and which showed positive niacin and nitrate reduction tests and a low catalase activity. Any growth seen within seven days of incubation or which did not meet the above criteria was considered as NTM [3].

The genotype MTBDR plus assay, a multiplex PCR DNA strip test assay (Hains Lifescience, Nehren, Germany) was done to detect the molecular genetic identification of *M. tuberculosis* complex from BA culture positive samples. Step by step procedure was done as instructed by the manufacturers (DNA extraction, multiplex amplification with biotinylated primers and reverse hybridisation) [9]. The study was approved by the Ethics Committee, Osmania Medical College, Koti, Hyderabad, India.

STATISTICAL ANALYSIS

OpenEpi: Open Source Epidemiologic Statistics for Public Health, Version. www.OpenEpi.com, updated 2013/04/06 was used for statistical analysis. The Yates corrected Chi-square test was used to analyse the given data. A p-value <0.05 was considered significant.

RESULTS

A total of 100 cases were included in the study, 70 were males and 30 females. The positivity rate of ZN stain was found to be maximum in the 20-40 age group (11 of 21 smear positives). Out of 100 cases, 14 were HIV positive of which five were smear positive for AFB. About 21 were smear positive for AFB of which 19 showed growth on BA, 16 on LJ medium and 21 on BACTEC MGIT 960. Out of the 79 smear negative cases, 15 showed growth on BA, 14 on LJ and 17 on MGIT 960. *M. tuberculosis* was isolated from 39

samples overall, sensitivities being 87.17% (34 out of 39) for BA, 76.9% (30 out of 39) for LJ and 97.3% (38 out of 39) for BACTEC MGIT 960.

Total No. of Samples (n=100)	BA (n=34)	LJ (n=30)	MGIT (n=38)	Any culture positive (n=39)
Origin of Samples				
Sputum (n=95)	32	29	36	37
Pleural fluid (n=3)	1	1	1	1
Pus (n=1)	1	0	1	1
Lymph node (n=1)	0	0	0	0
Smear positive (n=21)				
Number of culture positive	19	16	21	21
Contaminated samples	0	0	0	0
Non <i>Tuberculous Mycobacteria</i>	0	0	0	0
Smear Negative (n=79)				
Number of culture positive	15	14	17	18
Contaminated samples	2	6	12	-
Non <i>Tuberculous Mycobacteria</i>	3	2	2	-
Sensitivity	87.17%	76.9%	97.3%	-

[Table/Fig-1]: Performance of different culture types.

About 39 *M. tuberculosis* isolates were obtained from any of the 3 culture media. NTM isolates were 3 on BA, 2 on LJ and 2 on MGIT 960 and the rate of contamination was 12% in MGIT 960, 6% on LJ and 2% in BA [Table/Fig-1].

Variables	Blood agar when compared to Lowenstein Jensen medium	Blood agar when compared to BACTEC MGIT 960.
Specificity	100%	97.06%
Positive predictive value	100%	92.42%
Negative predictive value	94.29%	86.84%

[Table/Fig-2]: Specificity, positive predictive value, negative predictive value of Blood agar when compared to Lowenstein Jensen and BACTEC MGIT 960.

In our study, the specificity, positive predictive value, and negative predictive value of BA were 97.06%, 92.42% and 86.84% when compared to BACTEC MGIT 960. The specificity, positive predictive value, and negative predictive value of BA when compared to LJ was observed to be 100%, 100%, and 94.29%, respectively [Table/Fig-2].

Variables	Blood agar culture positive	Blood agar culture negative
LJ culture positive	30	0
LJ culture negative	4	66

[Table/Fig-3]: Comparison of results of culture between Blood agar and Lowenstein Jensen medium. $\chi^2 = 79.04$; p value = <0.001.

In the present study, a significant difference (p-value < 0.001), was observed to detect *M. tuberculosis* culture positivity on BA when compared to LJ [Table/Fig-3]. Thirty isolates were culture positive on both BA and LJ medium, in addition BA could isolate four cases which were negative on LJ medium [Table/Fig-3].


Variables	Blood agar culture positive	Blood agar culture negative
BACTEC MGIT 960 positive	33	5
BACTEC MGIT 960 negative	1	61

[Table/Fig-4]: Comparison of results of culture between Blood agar and BACTEC MGIT 960. $\chi^2 = 72.51$; p value = <0.001

In the present study, a significant difference (p -value <0.001) was observed to detect *M. tuberculosis* culture positivity on BA and MGIT 960. MGIT could isolate additional 5 *M. tuberculosis* which were culture negative on BA. But, there was one isolate which was negative by BACTEC MGIT 960 but positive for culture on BA [Table/Fig-4].

Duration of incubation (days)	Blood agar – No. of isolates	LJ – No. of isolates	BACTEC MGIT 960-No. of isolates
1 – 7	-	-	2
8 – 14	25	-	30
15 – 21	9	2	6
22 – 28	-	17	-
29 – 35	-	10	-
36 – 42	-	1	-
Total	34	30	38

[Table /Fig-5]: Duration of incubation for isolation of *Mycobacterium tuberculosis* on Blood Agar, Lowenstein Jensen and MGIT 960.



[Table/Fig-6]: Blood agar showing growth of *Mycobacterium tuberculosis*.

The minimum time taken by *M. tuberculosis* to grow was eight days on BA, 21 days on LJ and six days on MGIT 960. The maximum time taken was 18 days on BA, 36 days on LJ, and 16 days on MGIT 960 [Table/Fig-5].

Out of the 37 culture positive cases on BA, 34 samples showed TUB band in the Genotype MTBDR plus assay version 2.0, thus confirming the isolates belonging to MTB complex and remaining three showed no band confirming it as NTM.

DISCUSSION

Robert Koch isolated *M. tuberculosis* from freshly crushed pulmonary tubercles using heat coagulated sheep and beef serum medium after 10 days of incubation in tubes [10]. Isolates were recovered using BA, but were replaced by an egg based agar which became the recommended medium for isolation of *M. tuberculosis* in 1920's as the egg based medium is easy to sterilize [10]. Drancourt M et al., were the first to report the growth of *M. tuberculosis* colonies on BA [10]. The Indian studies conducted by Mathur ML et al., Rajdev S et al., Palange P et al., also have positively evaluated the utility of BA for primary isolation as well as Drug Susceptibility Testing (DST) of *M. tuberculosis* [4,11,12].

In our study, BA showed sensitivity of 87.17% as comparable to findings of Mathur ML et al., 94.2%, Satti L et al., 90%, Palange P et al., 89.3%, Shidiki A and Pokhrel N 98.9% [4,6,12,13]. LJ medium and MGIT 960 showed sensitivity of 76.9% and 97.3% in the present study. The mean duration of isolation of *M. tuberculosis* on BA was 13.2 days which are similar to findings done by Mathur ML et al., 13.6 days and, Satti L et al., 14 days [4,6], Palange P et al., 18.4 days and Shidiki A and Pokhrel N reported 13.6±5.2 days on BA [12,13]. The mean duration of isolation in LJ was 26 days and in MGIT 960 was 11.8 days in our study. The colonies on BA are pinpoint in size, light gray, glistening and easily recognized when examined under bright light, gradually increase in size, become irregular, warty, or cauliflower like colony [Table/Fig-6]. Growth of NTM also observed on BA [Table/Fig-7]. But the limitation is that Mycobacteria species are morphologically indistinguishable on BA.

The cost to prepare one BA bottle was Rs 11/-, while one LJ medium was Rs 9/- and MGIT cost Rs 190/- per bottle. Also, it is easy to prepare in comparison to LJ.

In our study, BA showed good results when compared to gold standard LJ medium in sensitivity being 87.17%, mean duration of isolation being 13.2 days and least contamination rate being 2% while on LJ medium, sensitivity was 76.9%, mean duration of isolation was 26 days and contamination rate was 6%. Apart from our study, there are other 3 studies which compared BA with automated system for isolation of *M. tuberculosis*. Drancourt M



[Table/Fig-7]: Growth of Non *Tuberculosis Mycobacteria* on Blood Agar.

et al., compared BA with BACTEC 9000 MB broth and reported superior performance of BA (98.9% versus 92.6%) over BACTEC 9000 MB broth and duration of isolation was 26 days on BACTEC and 19 days on BA [5]. Mathur ML et al., reported isolation of *M. tuberculosis* took 9 days by MGIT and 13.6 days on BA [4]. Coban AY et al., reported 10 days on MGIT and 14 days on BA [7]. In our study, the contamination rate on BA was 2%, preferably least when compared to LJ and MGIT 960 which were 6% and 12%. Mathur ML et al., observed a contamination rate of 1.6% on BA and 7.8% on LJ medium [4]. Satti L et al., reported contamination rates of 2.8% and 4.2% on BA and LJ medium and Palange et al., showed 7.2% on BA and 4.8% on LJ medium [6,12]. No contamination was observed in BA, LJ or MGIT 960 supporting NTM growth.

Few studies like Coban AY et al., [14] and Shidiki A and Pokhrel among others have also gone to a further extent of using BA for (DST) of *M. tuberculosis* [13]. Though our study has not been involved in DST of *M. tuberculosis* but it can suggest that BA can support the growth of *M. tuberculosis*, the contamination rate was lower and would be at least as good as, if not better than, LJ medium for recovering *M. tuberculosis*.

LIMITATION

In the present study, the sample size was low. Growth of *M. tuberculosis* is morphologically indistinguishable from other Mycobacterial species on BA.

CONCLUSION

Based on the above findings, it can be said that BA can provide as a good substitute over LJ medium in terms of rate of isolation, duration of isolation and sensitivity. The contamination rate on BA was also low as compared to LJ medium. Although the sample size was less and BA lacked a little behind over MGIT 960 in terms of sensitivity and duration of isolation but the contamination rate was higher in MGIT as compared to BA. Also MGIT 960 requires sophisticated equipment, expensive kits and skilled manpower while BA is easy to prepare, inexpensive and requires only an incubator, basic equipment in clinical laboratory provided appropriate microbiological safety measures are in place to counter the highly infectious nature of *M. tuberculosis*.

REFERENCES

- [1] Govt. of India. Ministry of Health and Family welfare. Revised National TB Control Programme-annual status report: TB India; 2017;page 173.
- [2] Sharma K, Appannanavar SB, Goyal K, Sharma A. Recent advances in the diagnosis of tuberculosis. J Postgrad Med Edu Res. 2013;47(4):181-87.
- [3] Revised National TB Control Programme Training Manual for *Mycobacterium tuberculosis* Culture and Drug susceptibility testing. Central TB Division, Directorate General of Health Services. Printed April 2009. Accessed: 25th March, 2014: <http://tbcindia.nic.in/WriteReadData/1892s/6995271860Training%20manual%20M%20tuberculosis%20C%20DST.pdf>.
- [4] Mathur ML, Gaur J, Sharma R, Solanki A. Rapid culture of *Mycobacterium tuberculosis* on blood agar in resource limited setting. Dan Med Bull. 2009;56:208-10.
- [5] Drancourt M, Raoult D. Cost-effectiveness of blood agar for isolation of mycobacteria. PLoS Negl Trop Dis. 2007;1(2):e83.
- [6] Satti L, Ikram A, Coban AY, Martin A. Rapid direct testing of susceptibility of *Mycobacterium tuberculosis* to Isoniazide and Rifampicin on Nutrient and Blood agar in Resource limited settings. J Clin Microbiol. 2012;50(5):1659-62.
- [7] Coban AY, Akqunes A, Durupinar B. Evaluation of Blood agar medium for the growth of mycobacteria. Mikrobiyol Bul. 2011;45(4):617-22.
- [8] Koneman EW, Allen SD, Janda WM, Schrenkenberger PC, Winn WC, Procop GW, et al. Color Atlas and Textbook of Diagnostic Microbiology. 6th ed. Philadelphia: Lippincott Williams and Wilkins; 2006. Chapter 19, Mycobacteria; pp.1066-1124.
- [9] Hains Life Sciences. Genotype MTBDR plus – Manual of Standard Operating Procedures.
- [10] Drancourt M, Carrieri P, Gévaudan MJ, Raoult D. Blood agar and *Mycobacterium tuberculosis*: the end of a dogma. J Clin Microbiol. 2003;41(4):1710-11.
- [11] Rajdev S, Patel A, Mulla S. Use of blood agar along with Lowenstein Jensen media for rapid isolation of *Mycobacterium tuberculosis* for early drug susceptibility testing. J Acad Clin Microbiol. 2014;16(2):100-103.
- [12] Palange P, Narang R, Kandi V. Evaluation of culture media for isolation of *mycobacterium* species from human clinical specimens. Cureus. 2016;8(8):e7575.
- [13] Shidiki A, Pokhrel N. A comparison of blood and egg based media for the rapid isolation and drug susceptibility testing of *Mycobacterium Tuberculosis*. Research Journal of Pharmaceutical, Biological and Chemical Sciences. 2013;4(1); 736-43: 0975-8585.
- [14] Coban AY. Blood agar validation for susceptibility testing of isoniazid, rifampicin, ethambutol, and streptomycin to *Mycobacterium tuberculosis* isolates. PLoS ONE. 2013;8(2):e55370. Available from <https://doi.org/10.1371/journal.pone.0055370>.

PARTICULARS OF CONTRIBUTORS:

1. Assistant Professor, Department of Microbiology, RVM Institute of Medical Sciences, Hyderabad, Telangana, India.
2. Professor and Head, Department of Microbiology, Osmania Medical College, Hyderabad, Telangana, India.
3. Assistant Professor, Department of Microbiology, Osmania Medical College, Hyderabad, Telangana, India.

NAME, ADDRESS, E-MAIL ID OF THE CORRESPONDING AUTHOR:

Dr. Swetha Kamble,
Flat No.304, Meghana Homes, Near D.A.V Public School, Santoshima Nagar, Saffilguda, Malkajgiri, Secunderabad-500056,
Hyderabad, Telangana, India.
E-mail: swethakamble@gmail.com

Date of Submission: **Apr 03, 2017**

Date of Peer Review: **May 09, 2017**

Date of Acceptance: **Aug 23, 2017**

Date of Publishing: **Jan 01, 2018**

FINANCIAL OR OTHER COMPETING INTERESTS: None.